

KATANA™ Zirconia



New

90 Minutes Sintering
STML Shade Availability Expansion


**MULTI
LAYERED**

The logo "MULTI LAYERED" expresses original production technology of KATANA™ multi-layered materials.

UTML Ultra Translucent Multi Layered

STML Super Translucent Multi Layered


ML Multi Layered **HT*** High Translucent

TOTAL SINTERING TIME: 90 MINUTES

KATANA™ Zirconia, the innovative multi-layered zirconia, known for specialized fabrication of Full Contour Zirconia prosthetics is now speed capable. The AUSTROMAT® 674i® fast sintering furnace can now be used with KATANA™ Zirconia in the speed sintering process*¹. You can successfully finish sintering KATANA™ Zirconia in approximately 90 minutes*².

*¹ Only for single unit *² The total sintering time may be affected depending on: furnace condition; surrounding temperature; voltage; and electrical current, etc.

Fast sintering furnace




AUSTROMAT® 674i® (DEKEMA Dental-Keramiköfen GmbH)

Sintering schedule

STML/UTML series: 1560°C(2840°F)

ML/HT series: 1515°C(2759°F)




EXPANDING THE SHADE AVAILABILITY OF STML

STML exhibits an excellent balance of translucency and desired mechanical properties. Now available in 9 expanded shades for completion of a wider range of esthetic restorations.



STML Super Translucent Multi Layered



New

Printed color can be slightly different from actual color.

YOUR CONTACT

Kuraray Europe GmbH
BU Medical Products
Philipp-Reis-Straße 4
65795 Hattersheim am Main

Phone DACH
Phone Rest of Europe
Fax
E-Mail
Website

+49 (0) 69-305 35 835
+49 (0) 69-305 85 980
+49 (0) 69-305 98 35 835
dental.eu@kuraray.com
www.kuraraynoritake.eu



Kuraray Noritake Dental Inc.
300 Higashiyama, Miyoshi-cho, Miyoshi, Aichi 470-0293, Japan

*"KATANA" is a trademark of NORITAKE CO., LIMITED.